



Only Incubator Manufacturer Makes It Possible.



■ 2 Point Calibration with Standard Gas
CO₂/O₂ sensors can be calibrated.

■ Equipped with IR sensor & Galvanic Sensor. Accurate analysis with no influence from atmospheric change.

■ Portable
Compact, small design with built-in rechargeable battery.

■ Low Cost
Durable, highly reliable analysis at low price.

■ Easy Operation
Analysis can be completed with one single operation.

Model	AGA-2004		
CO ₂ Sensor	IR Sensor	O ₂ Accuracy	±0.1%
CO ₂ Range	0~18.0%	Weight	4.5kg
CO ₂ Accuracy	±0.1%	Power	AC100V 1A 50/60Hz
O ₂ Sensor	Galvanic Sensor	Outer Dimension	W150 x D280 x H187 (mm)
O ₂ Range	0~22.0%	Battery Life	8 hours depending on use