

# ASTECC Gas line Filters


ASTECC-6: Six-month Filter

ASTECC-4: Four-month Filter



## Achieve Better results with Gas Purified by Astec

Astec Gas Line Filters use a combination of activated carbon granules, a HEPA filter and pre-filters to remove VOCs and trap microorganisms, aldehydes and particles. Unique Micro Carbon Technology provides a large activated carbon surface area, enabling uniform gas flow and extremely high purifying capacity.



## Two types of Durability

Astec Gas Line filters are available in two variants: Astec Filter 6 (six month longevity) and the longer lasting Astec Filter 4 (four month longevity).



Line Connector : AS-082085



Holder : AS-082085

## High Quality

Manufactured and tested in Clean Room Environment, according to cGMP under strict rules for IVF use.

## Pure Materials for IVF use

Capsule and all internal materials are safe to use in IVF laboratories.

## No Chemicals Inside

Chemicals may have a negative effect on the embryo development and we have entirely eliminated the use of chemicals in our assembly and production line

## Quality Control

Each filter is tested in a Clean Room Environment with an Automated Quality Testing Device, which eliminates human errors entirely.

## Packing

Filters are individually packed and delivered as a set of four pieces in a sealed cardboard box which includes a certificate of conformity and instruction for use.



[ Please contact us below ]

Email : [astec-japan@astec-bio.com](mailto:astec-japan@astec-bio.com) TEL : +81-92-935-5666  
Or your local ASTEC distributor

[www.astec-bio.com/global/](http://www.astec-bio.com/global/)